



INSTRUCTIONS FOR USE LATEX CONTROLS

EN

Rev. 2.1bis – 2023-12-06

REF SI 305.100-A; SI 305.102-A; SI 305.300-A; SI 305.302-A

6/30 TESTS KIT



INTENDED USE

Latex Controls are intended for in-vitro diagnostic professional use only in the quantitative control of the calibration status of ALIFAX ESR line analyzers. The control is performed on three samples with different known turbidity values analyzed for transmittance related to ESR values.

DEVICE DESCRIPTION

Latex Control is provided in ready-to-use plastic capped tubes filled with 3ml aqueous suspension of synthetic latex at different known concentrations:

- 13x75mm Greiner Tube: REF: SI 305.100-A (6 test); SI 305.300-A (30 test)
- 11,5x66mm Sarstedt Tube: REF: SI 305.102-A (6 test); SI 305.302-A (30 test)

Each triplet allows up to 6 checks. The traceability is assured by the lot number and the barcode carried on each tube. The identification of the tubes is by an identification label indicating the rack (or rotor) loading position number and the barcode including the standard turbidimetric value.

PRINCIPLE OF THE METHOD

Good laboratory practices require the use of stable reference materials to verify the accuracy and precision of testing equipment and procedures. The technology applied by Alifax's ESR analyzers is Quantitative Capillary Photometry, which allows to obtain the ESR result of the sample, expressed in mm/hour, as per guidelines and reference method. For laboratory quality control needs it is then important to have at disposal a control system reproducible and easy to handle.

The Control kit is based on the use of three samples with known turbidity values, on which the analyzer performs transmittance measurements related to ESR values. The results obtained should fit the expected ranges. Otherwise the calibration of the instrument shall be verified.

MATERIALS PROVIDED

Each kit is made up of 1 or 5 triplets containing each the following 3 test-tubes:

- No.1 Latex Control level 2 ("LATEX Test tube L 2")
- No.1 Latex Control level 3 ("LATEX Test tube L 3")
- No.1 Latex Control level 4 ("LATEX Test tube L 4")

MATERIALS REQUIRED BUT NOT PROVIDED

Washing tubes containing 3ml distilled water

COMPATIBILITY

ALIFAX ESR line analyzers upgraded with software (SW) for latex management

- SW 6.01A onward: TEST1: (REF: SI 195.210/THL; SI 195.220/BCL; SI 195.230/SDL; SI 195.240/YDL; SI 195.250/MDL; SI 195.260/XDL); MicroTEST1 (REF: SI 199.101-LC); ROLLER 20 (REF: SI R20-LC)
- SW 1.0.0 onward: TEST1 2.0: (REF: SI 195.210/THL)
- SW 1.00A onward: ROLLER 20 (REF: SI R20-PN)
- SW 4.01A onward: ROLLER 20 (REF: SI R20-MC)
- SW 4.03A onward: ROLLER 10 PLUS NEEDLE (REF: SI R10-PN)
- SW 1.00.08 onward: JO-PLUS (REF: SI 804.100)

WARNINGS AND PRECAUTIONS

- For in-vitro diagnostic use only
- Do not use if the packaging is damaged.
- Do not use if product shows aggregations

- Handle with caution, avoiding ingestion, inhalation, contact with eyes, skin and clothes.
- Do not reuse the product more than six times

BENEFITS AND LIMITATIONS

The product can be used only with ALIFAX ESR line analyzers.

REACTIVE INGREDIENTS

Polymer based on styrene and butadiene.

PROCEDURE

The daily use of the control kit is recommended.

All kit components must be brought to room temperature prior to performing the test.

Perform the washing procedure according to the instrument operative manual.

Take three control tubes (for 30 test kit, verify that the three control tubes belong to the same column in the packaging) *and perform the control procedure according to the instrument operative manual.*

After the analysis of the control standards, the results obtained will be *displayed/printed* as ESR values (mm/h) in increasing order.

Only for SI R20-PN, check conformity of the 2 values (SENS1 and SENS2) for each level and if necessary calculate and record the average.

Store the kit at 4-8°C immediately after use.

RESULTS

The three ESR resulting values should fit the ranges reported both on the *display/print* itself and on the outer label in the table "REFERENCE VALUES". If the obtained values are outside the expected ranges, the calibration of the instrument shall be verified. In this case, call the technical service to recalibrate the instrument.

QUALITY CONTROL

Quality Control can be done by processing repeatability tests on the same instrument and monitoring the results obtained with turbidimetric standards on which the percentage Standard Deviation (% RSD) for each value is calculated. The average value must be under 10%.

STORAGE AND STABILITY CONDITIONS

The kit is stable until the expiration date shown on the label when stored in the original packaging, between 4-25°C. Do not use after the expiration date. Do not freeze the product (freezing can cause irreversible aggregations) or expose it to direct light during storage. The tubes are stable for 6 weeks since the first opening if properly stored at 4-8°C in the original packaging. No controlled temperature transport required provided that the temperature doesn't decrease under 0°C.

DISPOSAL

Reagents and samples disposal is under user's responsibility according to their specific features and to local law requirements.

BIBLIOGRAPHY

- Piva E, Pajola R, Tempurin V, Plebani M. A new turbidimetric standard to improve the quality assurance of the erythrocyte sedimentation rate measurement. Clinical Biochemistry 40 (2007) 491-495
- Choong-Hwan Cha, Young Joo Cha, Chan-Jeoung Park, Hyun Kyung Kim, Eun-Jong Cha, Duck Hee Kim, Honghooon,

ESR_IFU_SI305X0X_LATEXCONT_2-1bis_EN



ALIFAX S.r.l.

Via Merano,30 – 33045 Nimis (UD) - Italy














VAT number IT04337640280

www.alifax.com

Jae-Seol Jung, Mi-Jung Kim, Seongsoo Jang, Hyun-Sook Chi, Dong Soon Lee. Evaluation of the TEST 1 erythrocyte sedimentation rate system and intra- and inter-laboratory quality control using new latex control materials. Clin Chem Lab Med 2010;48(7):1043–1048

- Lapić I, Piva E, Spolaore F, Tosato F., Pelloso M, Plebani M. Automated measurement of the erythrocyte sedimentation rate: method validation and comparison. Clin Chem Lab Med 2019; 57(9): 1364–1373

SYMBOLS GLOSSARY

	Catalogue Number		For in vitro diagnostic use
	Manufacturer		Temperature Limit
	Use-by date yyyy-mm-dd		Batch Code
	Consult electronic instructions for use		Caution
	Kit content		Contains sufficient for <n> tests
	Keep away from sunlight		For US Only: Caution: US Federal Law restricts this device to sale by or on the order of the licensed practitioner.
	Do not use if package is damaged		

NOTICE TO THE USER [REGULATION (EU) 2017/746]

Any **serious incident** that has occurred in relation to the device shall be reported to the manufacturer ([email: vigilance@alifax.com](mailto:vigilance@alifax.com)) and the competent authority of the Member State in which the user and/or the patient is established.

Text written in italic blue identifies an addition or modification to the previous version